REF  
20020

# EASE-IT™

## Composite Heater

en

da

sv

no

fi

de

fr

pl

it

es

nl

sk

### INSTRUCTIONS FOR USE

#### CONTENT:

1 x Heating Unit with electronic temperature control.  
1 x Power Supply (Friwo FW7401M/12) with mains cable.

#### INTENDED USE:

The EASE-IT™ Composite Heater is designed for warming up dental composite materials prior to use in order to improve the handling properties.

#### ADDITIONAL EASE-IT™ PARTS:

**Jar** for warming unit-dose tips as well as composites placed directly in the inner plate.

**Inner plate**, sterilizable.

**Thermoblock** for pre-filled composite syringes. Turned upside-down, it can be used for unit-dose tips instead.

#### INSTRUCTIONS:

1. Temperature can be set between 40 – 60°C (104 – 140°F) according to chosen material and individual preference. Default setting is at 50°C (122°F).
2. Connect the Power Supply with mains cable to the Heating Unit and plug it into a 100-240V outlet. Place the unit-dose tips or composite syringes in the jar or the Thermoblock.
3. When the LED indicator turns from orange to green, the EASE-IT™ has reached the selected temperature.

The jar and the Thermoblock are both heat-insulating and can be removed from the heating unit and placed at work site.

#### CAUTION:

To avoid cross-contamination, do not put used (contaminated) unit-dose tips back into EASE-IT™ for re-use.

en

#### CLEANING / DISINFECTING:

All parts can be wiped over with a cloth with disinfectant. Do not immerse any parts into liquid. Do not autoclave.

#### WARRANTY:

The device is covered by a 2-year guarantee on materials and construction. In case of malfunction, please return the unit for repair over your EASE-IT™ dealer.

#### TECHNICAL DATA:

Voltage:	12VDC
Power:	14VA
Temperature range:	
Operating:	10°C – 35°C
Storage:	-20°C – 60°C
Humidity:	10 – 95%
Altitude:	Operation in 0-3000 m. Storage and transportation: No restrictions.
Humidity:	Operation, storage and transportation: 10-95%.

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EMC Technical Information

Table 1 Electromagnetic emissions

The "EASE-IT™" is intended for use in the electromagnetic environment specified below. The customer or the user of the "EASE-IT™" should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR11 IEC 60601-1-2 EN55011	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11 IEC 60601-1-2 EN55011	Group B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 60601-1-2 IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 60601-1-2 IEC 61000-3-3	Complies	

Table 2 Electromagnetic immunity

The "EASE-IT™" is intended for use in the electromagnetic environment specified below. The customer or the user of the "EASE-IT™" should assure that it is used in such an environment.

Immunity test	IEC 60601-1 Test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 60601-1-2 IEC61000-4-2	±6KV contact ±8KV air ±15KV air	±6KV contact ±8KV air ±15KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 40%.
Electrical fast transient/burst IEC 60601-1-2 IEC61000-4-4	±2KV for power supply lines ±1KV for input/output lines	±2KV for power supply lines ±1KV for input/output lines	Mains power supply quality should be that of typical residential area.
Surge IEC 60601-1-2 IEC61000-4-5	±1KV differential mode	±1KV differential mode	
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 60601-1-2 IEC61000-4-11	<5% UT for 0.5 cycle 40%UT for 5 cycles 70%UT for 25 cycles <5%UT for 250 cycles	<5% UT for 0.5 cycle 40%UT for 5 cycles 70%UT for 25 cycles <5%UT for 250 cycles	Mains power supply quality should be that of typical residential area.
Power frequency (50-60Hz) magnetic field. IEC 60601-1-2 IEC61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be that of typical residential area.

Table 3

## Electromagnetic immunity

The "EASE-IT™" is intended for use in the electromagnetic environment specified below. The customer or the user of the "EASE-IT™" should assure that it is used in such an environment.

Immunity test	IEC 60601-1 Test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 60601-1-2 IEC61000-4-6	10Vrms 150KHz to 80 MHz	10Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$
Radiated RF	3V/can	3V/m	$d = 1,2\sqrt{P}$ 80 MHz to Hz to 800 MHz
IEC 60601-1-2 IEC 61000-4-3	80MHz to 2,5GHz	±1KV differential mode	$d = 2,3\sqrt{P}$ 80 MHz to Hz to 2,5 GHz where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m).
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> .
			Interference may occur in the vicinity of equipment marked with the following symbol: 

*NOTE 1* At 80 MHz and 800 MHz, the higher frequency range applies.

*NOTE 2* These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

## Recommended separation distances between portable and mobile RF communications equipment and the "EASE-IT™"

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter (W)	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1,2\sqrt{P}$	80 MHz to 800 MHz $d=1,2\sqrt{P}$	800 MHz to 2,5 GHz $d=2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23


For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts

(W) according to the transmitter manufacturer.

*NOTE 1* At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

*NOTE 2* These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Explanation of symbols used in EASE-IT™ labelling.

	CE mark.
	Interference may occur in the vicinity of equipment marked with this symbol.
	To be disposed of as electronic waste according to WEEE Directive 2012/19/EU.
	Serial number.
	Manufacturing company.
	Read the instructions.
	Catalogue number.
	Temperature limitation.
	Humidity limitation.

## RØNVIG Dental Mfg. A/S

Gl. Vejlevej 59 • DK-8721 Daugaard • Denmark

Tel.: +45 70 23 34 11 • Fax: +45 76 74 07 98 • E-mail: [contact@ronvig.com](mailto:contact@ronvig.com)

[www.ronvig.com](http://www.ronvig.com)