SEE-IT™

en da no sv fi fr it es de nl cs lt

REF 6060 | 6061 | 6064

INSTRUCTIONS FOR USE

DESCRIPTION: SEE-IT™ is designed to identify and distinguish between infected and affected carious dentine. SEE-IT™ works by penetrating and staining the denatured collagen in the infected dentine, which is a by-product of the decay process. Even though affected dentine may be softer than normal healthy, hard dentine, it should not necessarily be removed before restoring the cavity. Affected carious dentine is normally able to remineralize. As SEE-IT™ only stains the infected carious dentine, it can be a valuable help in ensuring that all infected carious dentine is removed, permitting to preserve as much as possible of the underlying affected (remineralizable) dentine. However, SEE-IT™ does not preclude conventional tactile methods for the identification of carious dentine.

REF	Description
6060	SEE-IT™ Caries-detector Blue, 10 ml bottle
6061	SEE-IT™ Caries-detector Red, 10 ml bottle
6064	SEE-IT™ Caries-detector Blue, 5 x 1,4 ml syringes

INTENDED PURPOSE: SEE-IT™ is designed to identify and to distinguish between infected and affected carious dentine.

USER PROFILE: To be used only by dental professionals.

INGREDIENTS: Propylene Glycol
Hydroxypropyl Methylcellulose
Propyl Paraben
Purified Water
Red dye (E129; C.I. No. 16035) or blue dye (E 133; C.I. No. 42090)

BEFORE USE: 1. Check that the primary container is intact before use.
2. Applies to SEE-IT™ syringe (REF 6064): Make sure the needle and thread is correctly assembled before use to ensure that the needle doesn't fall off.

INSTRUCTIONS: 1. Wash and dry the carious dentine after examining carious dentine using conventional tactile methods.
2. Dispense 1-2 drops of SEE-IT™ into a mixing well or a dappen dish.
3. Use a small stick, sponge or mini brush for the application. Cover the carious dentine completely and rub gently.
4. Wait **exactly 10 seconds**.
5. Rinse with water while using the suction device.
6. Infected carious will now appear clearly stained.

After removal of the first stained layer, repeat the procedure until all the infected carious dentine has been removed and no staining is to be found. SEE-IT™ stains infected carious dentine also in difficult-to-see locations (e.g. under overhanging enamel) and can furthermore be used for the detection of difficult-to-find root canals.

CAUTION:  In case SEE-IT™ is spilled on cloth or skin, immediately wash with water and soap. SEE-IT™ is non-toxic and does not contain any harmful chemicals or dyes, but if accidentally spilled in the eye or an open wound, rinse with water immediately after SEE-IT™ is spilled and seek medical assistance. If SEE-IT™ is spilled in the eye, contact lenses must be removed prior to rinsing immediately after spilling of SEE-IT™.

RECOMMENDATIONS: 1. In the absence of sufficient data, the use during pregnancy and lactation is not recommended.
2. It is not recommended to use SEE-IT™ for children under 3 years.
3. It is recommended to wear gloves when applying the SEE-IT™.
4. Use suction device while rinsing to avoid swallowing of SEE-IT™.

STORAGE:  The SEE-IT™ can be used to the expiration date indicated on the packaging.
The product must be stored at 10-40°C/50-104°F when not in use.
Do not store near extreme heat sources or in direct sunlight.

DISPOSAL: Dispose of the SEE-IT™ if the product or packaging shows signs of damage.
After use, dispose of the SEE-IT™ must be handled in accordance with local regulations and instructions.

SERIOUS INCIDENT REPORTING: In case of any serious incident occurred in relation to the SEE-IT™ product it must be reported to RØNVIG Dental Mfg. A/S and the competent authority of the country where the user/patient is established.

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© RØNVIG Dental Mfg. A/S. ALL RIGHTS RESERVED. IFU_SEE-IT_multisprog_version04