FOCUSTM

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INDHOLDSFORTEGNELSE

FOCUS™ XTEND kontrolenhed side 4-5

LED lampeenheder side 6-7

FOCUS™ 5000 K

FOCUS™ 6500 K

FOCUS™ Nano 6500 K

Tabeller side 8-14

Symbolforklaring side 15

FOCUS™ XTEND
kontrolenhet

FOCUS™ 5000 K
FOCUS™ 6500 K

FOCUS™ Nano 6500 K

FOCUS™ XTEND KONTROLENHED

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Table 1: Electromagnetic emissions

The FOCUS™ XTEND LED is intended for use in the electromagnetic environment specified below. The customer or the user of the "FOCUS™ XTEND LED" should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2: Electromagnetic immunity

The FOCUS™ XTEND LED is intended for use in the electromagnetic environment specified below. The customer or the user of the "FOCUS™ XTEND LED should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC61000-4-2	±6KV contact ±8KV air	±6KV contact ±8KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 40%
Electrical fast transient/burst IEC61000-4-4	±2KV for power supply lines ±1KV for input/output lines	±2KV for power supply lines ±1KV for input/output lines	Mains power supply quality should be that of typical residential area.
Surge IEC61000-4-5	±1KV differential mode ±2KV common mode.	±1KV differential mode ±2KV common mode	Power frequency magnetic fields should be that of typical residential area.
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC61000-4-11	5% UT for 0.5 cycle 40%UT for 5 cycles 70%UT for 25 cycles <5%UT for 5 seconds	5% UT for 0.5 cycle 40%UT for 5 cycles 70%UT for 25 cycles <5%UT for 5 seconds	Mains power supply quality should be that of typical residential area.
Power frequency (50- 60Hz) magnetic field. IEC61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of typical residential area.

Table 3: Electromagnetic immunity

The FOCUS™ XTEND LED is intended for use in the electromagnetic environment specified below. The customer or the user of the FOCUS™ XTEND LED should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC61000-4-6 Radiated RF IEC61000-4-3	3Vrms 150KHz to 80 MHz 3V/can 80MHz to 2,5GHz	3Vrms 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,17\sqrt{P}$ $d=1,17\sqrt{P}$ 80 MHz to Hz to 800 MHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1

At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a)

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b)

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4:
Recommended separation distances between portable and mobile RF communications equipment and the FOCUS™ XTEND LED

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output Rated maximum output	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,17\sqrt{P}$	80 MHz to 800 MHz $d = 1,17\sqrt{P}$	800 MHz to 2,5 GHz $d = 3,5\sqrt{P}$
0,01	0,12	0,12	0,35
0,1	0,37	0,37	1,11
1	1,2	1,2	3,5
10	3,7	3,7	11,1
100	12	12	35

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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RØNVIG Dental Mfg. A/S. ALL RIGHTS RESERVED. IFU_FOCUS_DA_Version05_2026