



# Instruments for Restoration



## COMPO-JECT™

Applicators for dental composite materials in manufacturer prefilled unit-dose tips for direct application

- Takes universal pre-filled unit-dose tips
- Easy extrusion, no matter viscosity
- Easy to load and unload
- The tip is safely fixed and can be rotated in all directions
- No piston drawback - just snap in and go!
- Easy to reach all areas of the oral cavity
- Made from high grade, acid resistant stainless steel
- Can be sterilized by any method
- 5-year guarantee on material and construction



The clip-in-lock permits quick and easy loading/unloading and secures the tip safely.



**RØNVIG**  
Dental Mfg. A/S

Gl. Vejlevej 59 · DK-8721 Daugaard · Denmark  
Tel.: (+45) 70 23 34 11 · Fax: (+45) 76 74 07 98  
e-mail: export@ronvig.com · www.ronvig.com