


Injection System for Dental Local Anaesthesia

# CALAJECT



[← Back](#)

[↓ Download](#)

The computer-assisted injection system Calaject facilitates painless local anaesthesia in dentistry. The anaesthetic solution is injected at low pressure and an adjusted flow rate, with a gentle acoustic signal indicating the actual flow. The touchscreen panel is particularly easy to operate. Various functions, such as aspiration, are controlled by a foot switch.

## Statement by the Jury

The successful combination of black and silver components underlines the elegance and high quality of Calaject.



reddot award 2014  
winner

## Credits

Manufacturer:  
RØNVIG Dental Mfg. A/S, Dagaard, Denmark

Design:  
3PART a/s, Denmark

[www.ronvig.com](http://www.ronvig.com)

[www.3part.dk](http://www.3part.dk)



reddot

[Product Design](#)  
[Brands & Communication](#)  
[About Red Dot](#)  
[Red Dot Network](#)

[Press](#)  
[Contact](#)  
[My Red Dot](#)  
[Imprint](#)

[Privacy\\_policy](#)



Subscribe to newsletter

E-Mail

Sign up

Copyright  
Red Dot 2021